




## AndraStent<sup>®</sup> L / XL / XXL



- High radial force Cobalt Chromium bare metal stent
- Developed for large vessels
- Large product range  
diameters from 6 up to 32 mm  
lengths from 13 up to 57 mm

# AndraStent®

- Material: Cobalt Chromium for a lower crimping profile with smaller introducer sheath and high radial force.
- Design: Semi-flexible Hybrid Cell Design for improved expansion of the Stent even in the biggest diameter
- Stent Recoil: < 2.8% for low oversizing of the vessels

## Order Information AndraStent® L / XL / XXL Order Codes

Order Code	Expansion Diameter Stent (mm)										Minimum Recommended Introducer Sheath
	6	8	12	15	18	20	22	25	28	32	
AS13L	13	12	11	10	8						7F*
AS17L	17	17	15	14	11						
AS21L	21	21	19	18	15						
AS26L	26	25	23	21	17						
AS13XL		13	12	12	11	10	9	8			8F*
AS17XL		17	16	16	15	14	13	11			
AS21XL		21	20	20	19	18	16	14			
AS26XL		26	24	23	22	21	19	17			
AS30XL		30	29	28	27	25	23	21			
AS35XL		35	33	31	29	27	25	24			
AS39XL		39	36	35	33	31	29	27			
AS43XL		43	40	38	36	34	32	29			
AS48XL		48	44	42	40	37	35	33			
AS57XL		56	52	50	47	44	42	40			
AS17XXL			17	16	16	15	15	14	13	10	9F*
AS21XXL			21	20	20	20	19	18	16	13	
AS26XXL			26	25	24	23	22	21	19	16	
AS30XXL			30	29	28	27	26	24	22	18	
AS35XXL			35	34	33	32	31	29	26	21	
AS39XXL			38	37	36	35	34	33	30	25	
AS43XXL			41	41	40	40	39	37	34	28	
AS48XXL			46	45	45	44	43	40	38	32	
AS57XXL			55	53	52	51	49	47	43	37	

\* The diameter of the introducer sheath depends on the balloon + 2F for the stent

